




Nestled 'N' Nursing

Right Angle Nipple Ruler

For Use During Virtual Pump Flange Sizing


Directions:

1. Cut out the Right Angle Nipple Ruler on the solid black lines
2. Align the ruler with the top and side of your nipple, allowing measurement of the nipple for width and height
3. Use camera to take picture to allow for virtually measuring your nipple to assist in **obtaining proper breast flange sizing**



Position the ruler as shown in the image at left. Be sure that the picture you take is taken straight on (not at an angle) to ensure proper measurement reading.

Need help?
Schedule a Virtual
Flange Fitting:
www.nestled-n-nursing.com



| -1 inch - |

Print this document on 8.5"
by 11" paper at 100% scale.
Check print sizing by
comparing the
measurement above.